

the Fungitell® Bulletin

volume 10, issue 4

Topic: (1→3)-β-D-GLUCAN AND COCCIDIOIDOMYCOSIS

Fungitell® Bulletins are intended as technical advisory communications and as such are disseminated to the general public in order to highlight the significance of (1→3)-β-D-Glucan on human health. These communications do not promote a specific drug, therapy nor make any representation or suggestion concerning the suitability or effectiveness of a particular drug or therapy in patients harboring (1→3)-β-D-Glucan. Fungitell® is an adjunct diagnostic assay to be utilized in conjunction with clinical signs and symptoms for the diagnosis of invasive fungal infection. Fungitell® is currently 510(k) cleared for the detection and quantification of (1→3)-β-D-Glucan in human serum and should be used and interpreted only in a manner consistent with the current Instructions for Use.

Discussion:

Coccidioidomycosis, commonly known as Valley Fever, presents a significant challenge to clinicians in endemic regions due to its varied and often non-specific symptoms, which can be easily misdiagnosed as other infectious and non-infectious lung diseases¹. Coccidioidomycosis is contracted through the inhalation of arthroconidia from the soil-dwelling fungi *Coccidioides immitis* and *Coccidioides posadasii*, affecting both immunocompromised and immunocompetent individuals in endemic areas such as the southwestern United States (California, Arizona), Central America, and South America. Recent Centers for Disease Control & Prevention (CDC) data show annual increases in coccidioidomycosis cases and expanding areas of endemicity². A pan-fungal biomarker such as (1→3)-β-D-glucan, may help clinicians determine whether a fungal infection is present in patients with symptoms of community-acquired pneumonia¹.

Coccidioidomycosis encompasses a spectrum of disease. Most cases are asymptomatic or present as mild respiratory infections or flu-like symptoms. Approximately 1% of cases progress to life-threatening extrapulmonary disseminated disease, typically in immunocompromised individuals^{3,4}. The gold standard for diagnosis includes culture and histopathology. Culture results typically take 4-5 days to yield positive results, but in some cases can take up to 2-3 weeks⁵. Serologic antibody testing is challenging because many assays are performed in reference laboratories, increasing turnaround time. Additionally, antibody-based assays have low sensitivity (21-56%) in patients who cannot mount a robust antibody response, leading to repeat testing or missed infections⁶.

Recently, Fungitell® has emerged as a promising adjunct diagnostic, offering a rapid assay to aid in the evaluation of patients suspected of coccidioidomycosis while awaiting *Coccidioides* antigen/antibody testing results⁷. Al-Obaidi *et al.* showed that combined *Coccidioides* serology and (1→3)-β-D-glucan testing in immunocompromised patients can aid in the initial diagnostic workup where coccidioidomycosis is highly suspected, with 82% of patients showing a positive result in at least one of the assays⁸.

Fungitell®

Bulletin Volume 10, issue 4
Publication Date: September 2024

Corporate Headquarters
Associates of Cape Cod, Inc.
124 Bernard E. Saint Jean Drive
East Falmouth, MA 02536 USA
Tel: (508) 540-3444
www.acciusa.com

United Kingdom/Europe
Associates of Cape Cod Int'l., Inc.
Unit 1 F/G/H Academy Business Park
Lees Road, Knowsley, Liverpool L33 7SA, UK
Tel: (44) 151-547-7444
www.acciuk.co.uk

(1→3)-β-D-GLUCAN AND COCCIDIOIDOMYCOSIS

Thompson *et al.* demonstrated that (1→3)-β-D-glucan testing is most beneficial for patients with severe coccidioidomycosis requiring hospitalization. These findings were supported by a subsequent small-scale study by Zangeneh *et al.*, which showed (1→3)-β-D-glucan positivity in over 90% of patients who were highly positive for *Coccidioides* antigen⁹. Although Thompson *et al.* reported a modest sensitivity of 44% in hospitalized coccidioidomycosis patients, this improved with disease severity: (1→3)-β-D-glucan positivity was 53% in acute cases and 83% in disseminated and coccidioidal meningitis cases⁷. Stevens *et al.* (2016) specifically evaluated sensitivity and specificity of Fungitell® in cerebrospinal fluid (CSF) samples in patients with suspected coccidioidal meningitis. In this study, (1→3)-β-D-glucan testing outperformed CSF antibody testing for diagnosing coccidioidal meningitis¹⁰.

While (1→3)-β-D-glucan values were higher in disseminated and coccidioidal meningitis cases than in acute infections, they did not correlate with serum or CSF coccidioidal IgG antibody titres^{7,10}. Notably, (1→3)-β-D-glucan was elevated in some hyperacute infection cases before antibody responses were detectable. For example, Thompson *et al.* found elevated (1→3)-β-D-glucan in three patient samples before positive antibody testing⁷, and Stevens *et al.* found elevated (1→3)-β-D-glucan in five CSF samples that were positive for coccidioidal antigen but had no detectable antibodies¹⁰, suggesting utility in early-stage infections.

In a single case report, Fungitell® was evaluated by Wilmes *et al.* for its potential to guide the duration of antifungal therapy. Although (1→3)-β-D-glucan levels generally correlated with disease resolution, further work is needed to evaluate the utility of (1→3)-β-D-glucan for directing optimal treatment duration in this group of patients¹¹.

By identifying fungal infections early, (1→3)-β-D-glucan testing can facilitate timely intervention in patients hospitalized with coccidioidomycosis, coccidioidal meningitis, or hyperacute coccidioidomycosis, where the infection progresses rapidly^{7,12}. The timing of (1→3)-β-D-glucan sampling in relation to the *Coccidioides* lifecycle may be important, as *Coccidioides* spherules contain roughly 60% (1→3)-β-D-glu-

can by dry weight, compared to arthroconidia, which contain only 20%. This lifecycle difference of (1→3)-β-D-glucan may explain the reduced sensitivity of Fungitell® with respect to coccidioidal infections compared to other invasive fungal infections^{7,13}.

In summary, in combination with traditional diagnostic methods for coccidioidomycosis, (1→3)-β-D-glucan testing with Fungitell® offers a rapid and useful adjunct, particularly in the context of severe disseminated and hyperacute cases where timely identification and treatment are critical.

References:

1. Poplin V, Smith C, Milsap D, Zabel L, Bahr NC. Diagnosis of Pulmonary Infections Due to Endemic Fungi. *Diagnostics*. 2021; 11(5):856. <https://doi.org/10.3390/diagnostics11050856>
2. Benedict K, McCotter OZ, Brady S, *et al.* Surveillance for Coccidioidomycosis — United States, 2011–2017. *MMWR Surveill Summ* 2019;68(No. SS-7):1–15. DOI: <http://dx.doi.org/10.15585/mmwr.mm6807a1>.
3. Galgiani JN, Ampel NM, Blair JE, *et al.* 2016 Infectious Diseases Society of America (IDSA) Clinical Practice Guideline for the Treatment of Coccidioidomycosis. *Clin Infect Dis*. 2016;63(6):e112–e146. doi:10.1093/cid/ciw360
4. Rush B, Zeigler J, Cheng MP, Hrymak C, Lother S. Outcomes with severe pulmonary coccidioidomycosis and respiratory failure in the United States. *J Assoc Med Microbiol Infect Dis Can*. 2023;8(1):40–48. Published 2023 Mar 1. doi:10.3138/jammi-2022-0028
5. Azar MM, Malo J, Hage CA. Endemic Fungi Presenting as Community-Acquired Pneumonia: A Review. *Semin Respir Crit Care Med*. 2020;41(4):522–537. doi:10.1055/s-0040-1702194
6. Blair JE, Coakley B, Santelli AC, Hentz JG, Wengenack NL. Serologic testing for symptomatic Coccidioidomycosis in immunocompetent and immunosuppressed hosts. *Mycopathologia* 2006; 162: 317–324.
7. Thompson GR 3rd, Bays DJ, Johnson SM, Cohen SH, Pappagianis D, Finkelman MA. Serum (1→3)-β-D-glucan measurement in coccidioidomycosis. *J Clin Microbiol*. 2012;50(9):3060–3062. doi:10.1128/JCM.00631-12
8. Al-Obaidi MM, Ayazi P, Shi A, Campanella M, Connick E, Zangeneh TT. The Utility of (1→3)-β-D-Glucan Testing in the Diagnosis of Coccidioidomycosis in Hospitalized Immunocompromised Patients. *J Fungi (Basel)*. 2022;8(8):768. Published 2022 Jul 25. doi:10.3390/jof8080768
9. Zangeneh TT, Malo J, Luraschi-Monjagatta C, *et al.* Positive (1-3) B-d-glucan and cross reactivity of fungal assays in coccidioidomycosis. *Med Mycol*. 2015;53(2):171–173. doi:10.1093/mmy/myu077
10. Stevens DA, Zhang Y, Finkelman MA, Pappagianis D, Clemons KV, Martinez M. Cerebrospinal Fluid (1,3)-Beta-d-Glucan Testing Is Useful in Diagnosis of Coccidioidal Meningitis. *J Clin Microbiol*. 2016;54(11):2707–2710. doi:10.1128/JCM.01224-16
11. Wilmes D, Schui D, Held J, *et al.* Disseminated coccidioidomycosis: Monitoring of serologic markers for treatment response. *Med Mycol Case Rep*. 2020;29:25–28. Published 2020 May 30. doi:10.1016/j.mmcr.2020.05.007
12. Stevens DA, Martinez M, Sass G, *et al.* Comparative Study of Newer and Established Methods of Diagnosing Coccidioidal Meningitis. *J Fungi (Basel)*. 2020;6(3):125. Published 2020 Aug 4. doi:10.3390/jof6030125
13. Viriyakosol S, Fierer J, Brown GD, Kirkland TN. Innate immunity to the pathogenic fungus *Coccidioides posadasii* is dependent on Toll-like receptor 2 and Dectin-1. *Infect Immun*. 2005;73(3):1553–1560. doi:10.1128/IAI.73.3.1553-1560.2005

Fungitell®

Bulletin Volume 10, issue 4
Publication Date: September 2024

Corporate Headquarters
Associates of Cape Cod, Inc.
124 Bernard E. Saint Jean Drive
East Falmouth, MA 02536 USA
Tel: (508) 540–3444
www.acciusa.com

United Kingdom/Europe
Associates of Cape Cod Int'l, Inc.
Unit 1 F/G/H Academy Business Park
Lees Road, Knowsley, Liverpool L33 7SA, UK
Tel: (44) 151–547–7444
www.acciuk.co.uk