

# the Fungitell® BULLETIN

volume 11, issue 1

Topic:

## Serum (1→3)-β-D-Glucan in critically ill patients and the intensive care unit

Martin Hoenigl, M.D.

Division of Infectious Diseases,  
ECMM Excellence Center for Medical Mycology,  
Department of Internal Medicine,  
Medical University of Graz, Graz, Austria.

### Discussion:

Invasive fungal infections, particularly those caused by *Candida*, *Aspergillus* and *Pneumocystis* species, pose a significant threat to critically ill patients, requiring prompt antifungal treatment. Serum 1,3-β-D-Glucan (BDG) testing, using the Fungitell® assay, has proven to be a reliable method for diagnosing candidemia and invasive candidiasis in ICU patients, patients with underlying hematologic malignancies [Prüller 2014] [Lamoth 2021], and other hospitalized patients [Ostrosky-Zeichner 2005]. One common characteristic feature shared by those fungi is the production of the polysaccharide BDG, a component of the fungal cell wall that serves as a diagnostic marker for invasive fungal infection. While BDG testing has a high negative predictive value (NPV), making it especially useful for safely discontinuing empirical antifungal therapy, it also plays a role in antifungal stewardship by helping reduce unnecessary treatment. This article explores the diagnostic role of BDG in managing invasive *Candida* infections in the ICU, its limitations with other fungi like *Aspergillus*, and its potential broader applications in gut health and emerging research areas.

*Fungitell® Bulletins are intended as technical advisory communications and as such are disseminated to the general public in order to highlight the significance of (1→3)-β-D-Glucan on human health. These communications do not promote a specific drug, therapy nor make any representation or suggestion concerning the suitability or effectiveness of a particular drug or therapy in patients harboring (1→3)-β-D-Glucan. Fungitell® is an adjunct diagnostic assay to be utilized in conjunction with clinical signs and symptoms for the diagnosis of invasive fungal infection. Fungitell® is currently 510(k) cleared for the detection and quantification of (1→3)-β-D-Glucan in human serum and should be used and interpreted only in a manner consistent with the current Instructions for Use.*

# Fungitell®

Publication Date: January 2025  
CORP\_0721

Corporate Headquarters  
Associates of Cape Cod, Inc.  
124 Bernard E. Saint Jean Drive  
East Falmouth, MA 02536 USA  
Tel: 508.540.3444  
www.acciusa.com

United Kingdom/Europe  
Associates of Cape Cod Int'l., Inc.  
Unit 1 F/G/H Academy Business Park  
Lees Road, Knowsley, Liverpool L33 7SA, UK  
Tel: (44) 151.547.7444  
www.acciuk.co.uk

**acc**   
PROTECTION  
THROUGH  
DETECTION



## Diagnosis in the ICU - *Candida*

Invasive *Candida* infections are life-threatening and require immediate initiation of antifungal treatment at the time of suspicion, to give patients a reasonable chance of survival and prevent serious complications such as *Candida* endocarditis [Clancy 2018] [Thompson 2023]. This observation naturally results in overtreatment, including empirical antifungal treatment of patients who in fact don't develop invasive *Candida* infections. The high sensitivity and negative predictive value (NPV) for invasive *Candida* infections, may make this marker an ideal candidate for utilization in antifungal stewardship efforts. Serum BDG has proven to be a reliable tool for antifungal treatment stratification [Posteraro 2016], particularly useful for guiding the cessation of empirical antifungal therapy [Sandhar 2022] [Bassetti 2024] due to its high NPV of ~90-95% [Bassetti 2024]. Given that serum BDG tests become negative under appropriate antifungal therapy (and are therefore also used for treatment monitoring) [Jaijakul 2012] it is important to obtain a sample for serum BDG testing before initiation of antifungal therapy, otherwise the NPV would decrease [Rautemaa-Richardson 2018]. Of note, sensitivity of serum BDG may be lower for diagnosing infections caused by *Candida auris* or *Candida parapsilosis* when compared to other *Candida* spp. [Mikulska 2016] [Farooqi 2021]. Multiple studies have demonstrated that negative serum BDG levels, if tested from samples drawn before initiation of systemic antifungal treatment, can be used to safely discontinue empirically initiated antifungal therapy in patients [Prattes 2014] [Nucci 2016].

Conversely, while serum BDG testing is valuable when used in a targeted approach, positive results from broad screening in ICU patients may not reliably justify initiating antifungal therapy in suspected cases of *Candida* infection. Data from a recent randomised controlled trials examining the use of antifungal therapy in ICU patients with sepsis at a high risk of *Candida* infection determined that, in this select group of at-risk patients, serum BDG guided initiation of antifungal therapy did not improve 28-day mortality and may be associated with overtreatment (i.e., treating patients with antifungals without invasive *Candida* infection) [Bloos 2022]. This suggests that the effectiveness of BDG as a diagnostic lies in its use within a more targeted approach. Serum BDG testing has also shown promise in diagnosing neonatal candidiasis [Cornu 2018] [Cohen 2020], and also *Pneumocystis jirovecii* pneumonia [Mercier 2019]. Of note, serum BDG has also proven reliable for predicting outcome in patients with candidemia [Agnelli 2020] [Esteves 2021]. The diagnostic potential of BDG in other patient materials remains limited, but CSF BDG testing has shown high sensitivity for diagnosing CNS *Candida* infections [Forster 2022], while high BDG titers in peritoneal fluid were associated with intra-abdominal candidiasis [Nourry 2023].

## Diagnosis in the ICU - *Aspergillus*

The role of serum BDG for diagnosing invasive aspergillosis in the ICU is limited. This is primarily due to the invasive growth pattern of *Aspergillus*, which tends to remain

localised in the airways, particularly in non-neutropenic ICU patients [Egger 2022], rather than disseminating into the bloodstream where BDG could be detected. As a result, serum BDG testing may not capture all cases of invasive pulmonary aspergillosis effectively. The challenges in diagnosing *Aspergillus* infections in non-neutropenic patients are well described [Colombo 2017], and while BDG may be useful as part of a broader strategy, its specificity limits utility for diagnosing *Aspergillus* infections. Furthermore, while BDG testing from urine has been explored as a potential adjunct to diagnosis, it has shown limited diagnostic potential [Raggam 2015], likely due to the variability in fungal dissemination and clearance patterns in the body. Consequently, alternative diagnostic tools, such as galactomannan combined with direct imaging, may be more effective for detecting invasive aspergillosis, especially in ICU patients.

## Additional Considerations

While the high sensitivity of BDG testing is considered the main strength in the ICU setting, a positive result only indicates fungal cell wall components in the blood, which can occur even without invasive *Candida* infection. Fungi, alongside bacteria, are major components of the gut microbiome [Hallen-Adams 2017], including *Candida* spp., *Saccharomyces* spp. and other yeasts, as well as molds like *Aspergillus* spp. and *Fusarium* spp. [Witherden 2017]. The role of the fungal gut microbiome, the mycobiome, has garnered increasing attention in recent years. Researchers are now focusing on the mycobiome's role and dynamics in health and disease, as well as the important interactions between its microorganisms - bacteria, viruses, archaea, fungi [Zhang 2022] [Egger 2024]. For example, Kusakabe and colleagues recently found elevated anti-*Candida* IgG antibodies in severe COVID-19 patients with intestinal dysbiosis, *Candida* overgrowth and systemic neutrophilia [Kusakabe 2023], indicating that the mycobiome may contribute to lung hyperinflammation. This finding highlights that gut dysbiosis of the mycobiome may be a contributing factor to the immunopathology of severe COVID-19 disease. Although fewer studies have evaluated the lung mycobiome's role in driving hyperinflammation, both *Candida* colonization of the lung and higher levels of BDG in bronchioalveolar lavage fluid, where BDG has limited diagnostic utility due to frequent *Candida* and *Pneumocystis* colonization [Mutschlechner 2015], have been shown to predict mortality in ICU patients in the absence of any association with lung infection [Reischies 2016]. The mechanisms underlying these associations remains unclear. However, in another context, the fungal cell wall component BDG has not only been shown to be elevated in serum from patients living with chronic HIV [Hoenigl 2019] [Farhour 2018], advanced liver cirrhosis [Egger 2023], and other conditions associated with microbial translocation, including after abdominal surgery [Szyszkowitz 2018], but has also been identified as a marker of immune activation and thereby a driver of persistent inflammation, that correlates with clinical outcomes [Hoenigl 2019] [Weiner 2019] [Yang 2017] [Gianella 2019].

## In Closing

In conclusion, serum BDG testing has a variety of important roles in the ICU population. Besides its diagnostic potential for invasive *Candida* and also *Pneumocystis* infections, serum BDG testing may be an effective tool for antifungal stewardship informing discontinuation of empirical antifungal therapy in patients that don't actually have *Candida* disease. For all indications a fast turnaround of the test, eg by utilizing the Fungitell STAT® assay [White 2021] or an automated Fungitell® assay protocol [Prüller 2014], can assure early discontinuation of antifungal therapy, thereby benefitting the patients and saving costs.

## Hoenigl Bio:

Martin Hoenigl, M.D., Assoc.

Prof., FECMM is an Associate Professor of Translational Mycology at the Division of Infectious Diseases, Medical University of Graz, Austria.

He has obtained his *venia docendi* in internal medicine in 2012, and is author to over 300 pub med listed publications in the field of infectious diseases, the majority in leading authorships (i.e. first or last author; ORCID: 0000-0002-1653-2824). Dr. Hoenigl has particular expertise in conducting research on clinical mycology, including fungal diagnostics and pharmacology of antifungal drugs and correlation with clinical findings. He is a past president of the European Confederation of Medical Mycology (ECMM), and the delegate of the Austrian Society for Medical Mycology (OEGMM). Dr. Hoenigl is also the founder of the ECMM Academy (together with Prof. Cornely) and serves as an associate editor for the journal *Mycoses* and deputy editor of *Mycopathologia*. Dr. Hoenigl has been awarded with the Researcher of the year 2011 award at the Medical University of Graz, and with the Research Promotion award 2014 of the German Speaking mycological society.



## References:

- Prüller F, Wagner J, Raggam RB, *et al.* Automation of serum (1-3)-beta-D-glucan testing allows reliable and rapid discrimination of patients with and without candidemia. *Med Mycol.* 2014;in press
- Lamoth F, Akan H, Andes D, *et al.* Assessment of the Role of 1,3-β-D-Glucan Testing for the Diagnosis of Invasive Fungal Infections in Adults. *Clinical infectious diseases: an official publication of the Infectious Diseases Society of America.* Mar 12 2021;72(Suppl 2):S102-s108. doi:10.1093/cid/ciaa1943
- Ostrosky-Zeichner L, Alexander BD, Kett DH, *et al.* Multicenter clinical evaluation of the (1-->3) beta-D-glucan assay as an aid to diagnosis of fungal infections in humans. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America.* 2005;41(5):654-659. doi:10.1086/432470
- Mikulska M, Giacobbe DR, Furfaro E, *et al.* Lower sensitivity of serum (1,3)-β-D-glucan for the diagnosis of candidaemia due to *Candida parapsilosis*. *Clinical microbiology and infection : the official publication of the European Society of Clinical Microbiology and Infectious Diseases.* Jul 2016;22(7):646.e5-8. doi:10.1016/j.cmi.2016.05.020
- Farooqi J, Niamatullah H, Irfan S, Zafar A, Malik F, Jabeen K. Comparison of β-D-Glucan levels between *Candida auris* and other *Candida* species at the time of candidaemia: a retrospective study. *Clinical microbiology and infection : the official publication of the European Society of Clinical Microbiology and Infectious Diseases.* Oct 2021;27(10):1519.e1-1519.e5. doi:10.1016/j.cmi.2021.05.031
- Bloos F, Held J, Kluge S, *et al.* (1-->3)-β-D-Glucan-guided antifungal therapy in adults with sepsis: the CandiSep randomized clinical trial. *Intensive Care Med.* Jul 2022;48(7):865-875. doi:10.1007/s00134-022-06733-x
- Cornu M, Goudjil S, Kongolo G, *et al.* Evaluation of the (1,3)-β-D-glucan assay for the diagnosis of neonatal invasive yeast infections. *Medical mycology.* Jan 1 2018;56(1):78-87. doi:10.1093/mmy/myx021
- Cohen JF, Ouziel A, Matczak S, *et al.* Diagnostic accuracy of serum (1,3)-beta-d-glucan for neonatal invasive candidiasis: systematic review and meta-analysis. *Clinical microbiology and infection : the official publication of the European Society of Clinical Microbiology and Infectious Diseases.* Mar 2020;26(3):291-298. doi:10.1016/j.cmi.2019.09.010
- Mercier T, Guldentops E, Patteet S, Beuselinck K, Lagrou K, Maertens J. Beta-d-Glucan for Diagnosing *Pneumocystis Pneumonia*: a Direct Comparison between the Wako b-Glucan Assay and the Fungitell Assay. *Journal of Clinical Microbiology.* 2019/05/24/ 2019;57(6)doi:10.1128/JCM.00322-19
- Agnelli C, Bouza E, Del Carmen Martínez-Jiménez M, *et al.* Clinical Relevance and Prognostic Value of Persistently Negative (1,3)-β-D-Glucan in Adults With Candidemia: A 5-year Experience in a Tertiary Hospital. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America.* Apr 15 2020;70(9):1925-1932. doi:10.1093/cid/ciz555
- Esteves P, Lopes Lima S, Salles de Azevedo Melo A, Maria Beirão E, Nucci M, Colombo AL. (1,3)-β-D-glucan is able to predict therapeutic failure of patients with candidemia and not only mortality. *Mycoses.* Mar 2021;64(3):264-271. doi:10.1111/myc.13224
- Egger M, Prüller F, Krause R, Prattes J, Hoenigl M. Utility of Serum 1,3-β-D-Glucan Testing for Diagnosis and Prognostication in COVID-19-Associated Pulmonary Aspergillosis. *Microbiol Spectr.* May 31 2022:e0137322. doi:10.1128/spectrum.01373-22
- Forster J, Hoenigl M, Suerbaum S, Wagener J, Dichtl K. Serologic biomarkers in *Candida* and *Aspergillus* infections of the central nervous system: a comparison of galactomannan, mannan, and β-1,3-D-gucan testing from serum and cerebrospinal fluid. *Mycoses.* May 4 2022;doi:10.1111/myc.13451
- Nourry É, Wallet F, Darien M, *et al.* Use of 1,3-Beta-D-Glucan concentration in peritoneal fluid for the diagnosis of intra-abdominal Candidiasis in Critically-ill patients. *Med Mycol.* Mar 20 2023;doi:10.1093/mmy/myad029
- Raggam RB, Fischbach LM, Prattes J, *et al.* Detection of (1-->3)-beta-D-glucan in same-day urine and serum samples obtained from patients with haematological malignancies. *Mycoses.* 2015;58(7):394-398. doi:10.1111/myc.12328 [doi]
- Hallen-Adams HE, Suhr MJ. Fungi in the healthy human gastrointestinal tract. *Virulence.* Apr 3 2017;8(3):352-358. doi:10.1080/21505594.2016.1247140
- Witherden EA, Shoaie S, Hall RA, Moyes DL. The Human Mucosal Mycobiome and Fungal Community Interactions. *J Fungi (Basel).* Oct 7 2017;3(4)doi:10.3390/jof3040056
- Zhang F, Aschenbrenner D, Yoo JY, Zuo T. The gut mycobiome in health, disease, and clinical applications in association with the gut bacterial microbiome assembly. *Lancet Microbe.* Sep 28 2022;doi:10.1016/S2666-5247(22)00203-8
- Egger M, Hoenigl M. One for All and All for One - Multi-Kingdom Interplay in Severe Viral Pneumonia. *Am J Respir Crit Care Med.* Jul 30 2024;doi:10.1164/rccm.202406-1237ED
- Kusakabe T, Lin WY, Cheong JG, *et al.* Fungal microbiota sustains lasting immune activation of neutrophils and their progenitors in severe COVID-19. *Nat Immunol.* Nov 2023;24(11):1879-1889. doi:10.1038/s41590-023-01637-4
- Mutschlechner W, Risslegger B, Willinger B, *et al.* Bronchoalveolar Lavage Fluid (1,3)beta-D-Glucan for the Diagnosis of Invasive Fungal Infections in Solid Organ Transplantation: A Prospective Multicenter Study. *Transplantation.* 2015;99(9):e140-4. doi:10.1097/TP.0000000000000635 [doi]



Reischies FM, Prattes J, Pruller F, *et al.* Prognostic potential of 1,3-beta-d-glucan levels in bronchoalveolar lavage fluid samples. *The Journal of infection.* 2016;72(1):29-35. doi:10.1016/j.jinf.2015.09.016 [doi]

Hoenigl M, Moser CB, Funderburg N, *et al.* Soluble Urokinase Plasminogen Activator Receptor Is Predictive of Non-AIDS Events During Antiretroviral Therapy-mediated Viral Suppression. *Clinical infectious diseases: an official publication of the Infectious Diseases Society of America.* 2019;69(4):676-686. doi:10.1093/cid/ciy966

Weiner L, Retuerto M, Hager C, *et al.* Fungal Translocation is Associated with Immune Activation and Systemic Inflammation in Treated HIV. *AIDS Res Hum Retroviruses.* Feb 20 2019;doi:10.1089/aid.2018.0252

Farhour Z, Mehraj V, Chen J, Ramendra R, Lu H, Routy JP. Use of (1-->3)-beta-d-glucan for diagnosis and management of invasive mycoses in HIV-infected patients. *Mycoses.* 2018;doi:10.1111/myc.12797 [doi]

Egger M, Horvath A, Prüller F, *et al.* Fungal translocation measured by serum 1,3-β-D-glucan correlates with severity and outcome of liver cirrhosis-A pilot study. *Liver Int.* Jun 19 2023;doi:10.1111/liv.15648

Szyszkowitz A, Zurl C, Herzeg A, *et al.* Serum 1,3-Beta-D-Glucan Values During and After Laparoscopic and Open Intestinal Surgery. *Open Forum Infect Dis.* Dec 2018;5(12):ofy296. doi:10.1093/ofid/ofy296

Yang AM, Inamine T, Hochrath K, *et al.* Intestinal fungi contribute to development of alcoholic liver disease. *The Journal of clinical investigation.* 2017;doi:10.1172/JCI90562 [doi]

Gianella S, Letendre SL, Iudicello J, *et al.* Plasma (1-->3)-beta-D-glucan and suPAR levels correlate with neurocognitive performance in people living with HIV on antiretroviral therapy: a CHARTER analysis. *J Neurovirol.* Jul 11 2019;doi:10.1007/s13365-019-00775-6

Clancy CJ, Pappas PG, Vazquez J, *et al.* Detecting Infections Rapidly and Easily for Candidemia Trial, Part 2 (DIRECT2): A Prospective, Multicenter Study of the T2Candida Panel. *Clin Infect Dis.* May 17 2018;66(11):1678-1686. doi:10.1093/cid/cix1095

Thompson GR, 3rd, Jenks JD, Baddley JW, *et al.* Fungal Endocarditis: Pathophysiology, Epidemiology, Clinical Presentation, Diagnosis, and Management. *Clin Microbiol Rev.* Sep 21 2023;36(3):e0001923. doi:10.1128/cmr.00019-23

Posteraro B, Tumbarello M, De Pascale G, *et al.* (1,3)-β-D-Glucan-based antifungal treatment in critically ill adults at high risk of candidaemia: an observational study. *J Antimicrob Chemother.* Aug 2016;71(8):2262-9. doi:10.1093/jac/dkw112

Sandhar TK, Chhina DK, Gupta V, Chaudhary J. Role of (1-3)-β-D-Glucan Test in the Diagnosis of Invasive Fungal Infections among High-Risk Patients in a Tertiary Care Hospital. *J Lab Physicians.* Sep 2022;14(3):306-311. doi:10.1055/s-0042-1742632

Bassetti M, Giacobbe DR, Agvald-Ohman C, *et al.* Invasive Fungal Diseases in Adult Patients in Intensive Care Unit (FUNDICU): 2024 consensus definitions from ESGCIP, EFISG, ESICM, ECMM, MSGERC, ISAC, and ISHAM. *Intensive Care Med.* Mar 21 2024;doi:10.1007/s00134-024-07341-7

Jaijakul S, Vazquez JA, Swanson RN, Ostrosky-Zeichner L. (1,3)-beta-D-glucan as a prognostic marker of treatment response in invasive candidiasis. *Clinical infectious diseases: an official publication of the Infectious Diseases Society of America.* 2012;55(4): 521-526. doi:10.1093/cid/cis456 [doi]

Rautemaa-Richardson R, Rautemaa V, Al-Wathiqi F, *et al.* Impact of a diagnostics-driven antifungal stewardship programme in a UK tertiary referral teaching hospital. *J Antimicrob Chemother.* Dec 1 2018;73(12):3488-3495. doi:10.1093/jac/dky360

Prattes J, Hoenigl M, Rabensteiner J, *et al.* Serum 1,3-beta-d-glucan for antifungal treatment stratification at the intensive care unit and the influence of surgery. *Mycoses.* 2014;57(11):679-686. doi:10.1111/myc.12221 [doi]

Nucci M, Nouér SA, Esteves P, *et al.* Discontinuation of empirical antifungal therapy in ICU patients using 1,3-β-D-glucan. *J Antimicrob Chemother.* Sep 2016;71(9):2628-33. doi:10.1093/jac/dkw188

White PL, Price JS, Backx M. Evaluation of the Performance of the Associates of Cape Cod, Inc. STAT Assay for the Diagnosis of Invasive Fungal Disease in Critical-Care Patients with COVID-19. *J Clin Microbiol.* Aug 18 2021;59(9):e0086921. doi:10.1128/jcm.00869-21

Colombo AL, de Almeida Júnior JN, Slavin MA, Chen SC, Sorrell TC. Candida and invasive mould diseases in non-neutropenic critically ill patients and patients with haematological cancer. *Lancet Infect Dis.* 2017;17(11):e344-e356. doi:10.1016/S1473-3099(17)30304-3